

Product Datasheet

Mouse Glypican 4 ELISA Kit (orb1212466)

Description	Mouse Glypican 4, GPC4 ELISA Kit
Reactivity	Mouse
Range	20-7000ng/L
Tested Applications	ELISA
Target	GPC4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	62,586
Uniprot ID	P51655
Solubility (25°C)	8000ng/L
Sensitivity	8.98ng/L
Entrez	14735
Expiration Date	Please enquire.