

Product Datasheet

Mouse Neuropilin-1 ELISA Kit (orb1212469)

Description	Mouse Neuropilin-1, NRP1 ELISA Kit
Reactivity	Mouse
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	NRP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	103,000
Uniprot ID	P97333
Solubility (25°C)	16ng/ml
Sensitivity	0.023ng/ml
Entrez	18186
Expiration Date	Please enquire.