

Product Datasheet

Mouse Elastin ELISA Kit (orb1212474)

Description	Mouse Elastin, ELN ELISA Kit
Reactivity	Mouse
Range	1-400ng/mL
Tested Applications	ELISA
Target	ELN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	71,938
Uniprot ID	P54320
Solubility (25°C)	480ng/mL
Sensitivity	0.57ng/mL
Entrez	13717
Expiration Date	Please enquire.