

Product Datasheet

Mouse Cartilage Oligomeric Matrix Protein ELISA Kit (orb1212479)

Description	Mouse Cartilage Oligomeric Matrix Protein, COMP ELISA Kit
Reactivity	Mouse
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	COMP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	82,342
Uniprot ID	Q9R0G6
Solubility (25°C)	240ng/ml
Sensitivity	0.22ng/ml
Entrez	12845
Expiration Date	Please enquire.