

Product Datasheet

Mouse Agouti Related Protein ELISA Kit (orb1212480)

Description	Mouse Agouti Related Protein, AGRP ELISA Kit
Reactivity	Mouse
Range	7-1500ng/L
Tested Applications	ELISA
Target	AGRP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,432
Uniprot ID	P56473
Solubility (25°C)	1600ng/L
Sensitivity	3.53ng/L
Entrez	11604
Expiration Date	Please enquire.