

Product Datasheet

Mouse Galectin 9 ELISA Kit (orb1212487)

Description	Mouse Galectin 9, GAL9 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	LGALS9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,036
Uniprot ID	008573
Solubility (25°C)	2400ng/L
Sensitivity	4.86ng/L
Entrez	16859
Expiration Date	Please enquire.