

Product Datasheet

Mouse Galectin 8 ELISA Kit (orb1212488)

Description	Mouse Galectin 8, GAL8 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	LGALS8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,162
Uniprot ID	Q9JL15
Solubility (25°C)	2400ng/L
Sensitivity	5.56ng/L
Entrez	56048
Expiration Date	Please enquire.