

Product Datasheet

Mouse Cathepsin K ELISA Kit (orb1212490)

Description	Mouse Cathepsin K, CTSK ELISA Kit
Reactivity	Mouse
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	CTSK
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,889
Uniprot ID	P55097
Solubility (25°C)	96ng/mL
Sensitivity	0.153ng/mL
Entrez	13038
Expiration Date	Please enquire.