

Product Datasheet

Mouse Very Low Density Lipoprotein Receptor ELISA Kit (orb1212491)

Description	Mouse Very Low Density Lipoprotein Receptor, VLDLR ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	VLDLR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	96,373
Uniprot ID	P98156
Solubility (25°C)	48ng/mL
Sensitivity	0.056ng/mL
Entrez	22359
Expiration Date	Please enquire.