

Product Datasheet

Mouse Fibroblast Growth Factor 15 ELISA Kit (orb1212507)

Description	Mouse Fibroblast Growth Factor 15, FGF15 ELISA Kit
Reactivity	Mouse
Range	2-600ng/L
Tested Applications	ELISA
Target	FGF15
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,236
Uniprot ID	035622
Solubility (25°C)	640ng/L
Sensitivity	1.01ng/L
Entrez	14170
Expiration Date	Please enquire.