

Product Datasheet

Mouse Caspase-9 ELISA Kit (orb1212528)

Description	Mouse Caspase-9 ELISA Kit
Reactivity	Mouse
Range	7-1500ng/L
Tested Applications	ELISA
Target	CASP9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,979
Uniprot ID	Q8C3Q9
Solubility (25°C)	1600ng/L
Sensitivity	3.29ng/L
Entrez	12371
Expiration Date	Please enquire.