

Product Datasheet

Mouse High Mobility Group Protein B1 ELISA Kit (orb1212540)

Description	Mouse High Mobility Group Protein B1, HMGB-1 ELISA Kit
Reactivity	Mouse
Range	1-400ng/ml
Tested Applications	ELISA
Target	HMGB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,894
Uniprot ID	P63158
Solubility (25°C)	480ng/ml
Sensitivity	0.55ng/ml
Entrez	15289
Expiration Date	Please enquire.