

Product Datasheet

Mouse Aquaporin 4 ELISA Kit (orb1212541)

Description	Mouse Aquaporin 4, AQP-4 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	AQP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,436
Uniprot ID	P55088
Solubility (25°C)	24ng/mL
Sensitivity	0.027ng/mL
Entrez	11829
Expiration Date	Please enquire.