

Product Datasheet

Mouse Granzymes B ELISA Kit (orb1212544)

Description	Mouse Granzymes B, GZMS-B ELISA Kit
Reactivity	Mouse
Range	7-1500ng/L
Tested Applications	ELISA
Target	GZMB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,470
Uniprot ID	P04187
Solubility (25°C)	1600ng/L
Sensitivity	3.26ng/L
Entrez	14939
Expiration Date	Please enquire.