

Product Datasheet

Mouse Ciliary Neurotrophic Factor ELISA Kit (orb1212556)

Description	Mouse Ciliary Neurotrophic Factor, CNTF ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	CNTF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,587
Uniprot ID	P51642
Solubility (25°C)	960ng/L
Sensitivity	1.62ng/L
Entrez	12803
Expiration Date	Please enquire.