

Product Datasheet

Mouse Transforming Growth Factor A ELISA Kit (orb1212573)

Description	Mouse Transforming Growth Factor A, TGF-A ELISA Kit
Reactivity	Mouse
Range	1-100ng/L
Tested Applications	ELISA
Target	TGFA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17018
Uniprot ID	P48030
Solubility (25°C)	128ng/L
Sensitivity	0.64ng/L
Entrez	21802
Expiration Date	Please enquire.