

Product Datasheet

Mouse Catalase ELISA Kit (orb1212594)

Description	Mouse Catalase, CAT ELISA Kit
Reactivity	Mouse
Range	0.2-70ng/ml
Tested Applications	ELISA
Target	CAT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,795
Uniprot ID	P24270
Solubility (25°C)	80ng/ml
Sensitivity	0.093ng/ml
Entrez	12359
Expiration Date	Please enquire.