

Product Datasheet

Mouse Insulin Like Growth Factor Binding Protein 4 ELISA Kit (orb1212609)

Description	Mouse Insulin Like Growth Factor Binding Protein 4, IGFBP4 ELISA Kit
Reactivity	Mouse
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	IGFBP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,807
Uniprot ID	P47879
Solubility (25°C)	96ng/mL
Sensitivity	0.154ng/mL
Entrez	16010
Expiration Date	Please enquire.