

Product Datasheet

Human Phosphatidylinositol 3,4 ELISA Kit (orb1212640)

Description	Human Phosphatidylinositol 3,4, 5- Trisphosphate-dependent Rac Exchanger 1 Protein, PREX1 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	PREX1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	186,203
Uniprot ID	Q8TCU6
Solubility (25°C)	32ng/ml
Sensitivity	0.027ng/ml
Entrez	57580
Expiration Date	Please enquire.