

Product Datasheet

Human Tumor Necrosis Factor Receptor Superfamily Member 13B ELISA Kit (orb1212652)

Description	Human Tumor Necrosis Factor Receptor Superfamily Member 13B, TNFRSF13B ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	TNFRSF13B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	31,816
Uniprot ID	O14836
Solubility (25°C)	1600ng/L
Sensitivity	3.32ng/L
Entrez	23495
Expiration Date	Please enquire.