

Product Datasheet

Human Hormone Sensitive Lipase ELISA Kit (orb1212656)

Description	Human Hormone Sensitive Lipase, HSL ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	LIPE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	116,598
Uniprot ID	Q05469
Solubility (25°C)	96ng/ml
Sensitivity	0.18ng/ml
Entrez	3991
Expiration Date	Please enquire.