

Product Datasheet

Human Lysophosphatidylcholine Acyltransferase 1 ELISA Kit (orb1212668)

Description	Human Lysophosphatidylcholine Acyltransferase 1, LPCAT1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	LPCAT1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,151
Uniprot ID	Q8NF37
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	79888
Expiration Date	Please enquire.