

Product Datasheet

Human Membrane Primary Amineoxidase ELISA Kit (orb1212677)

Description	Human Membrane Primary Amineoxidase, AOC3 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	AOC3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	84,622
Uniprot ID	Q16853
Solubility (25°C)	2400ng/L
Sensitivity	6.02ng/L
Entrez	8639
Expiration Date	Please enquire.