

Product Datasheet

Human Natural Killer Cell Receptor 2B4 ELISA Kit (orb1212679)

Description	Human Natural Killer Cell Receptor 2B4, CD244 ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	CD244
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,616
Uniprot ID	Q9BZW8
Solubility (25°C)	96ng/ml
Sensitivity	0.13ng/ml
Entrez	51744
Expiration Date	Please enquire.