

Product Datasheet

Human Aryl Hydrocarbon Receptor Nuclear Translocator-like Protein 1 ELISA Kit (orb1212694)

Description	Human Aryl Hydrocarbon Receptor Nuclear Translocator-like Protein 1, BMAL1 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	ARNTL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	68,762
Uniprot ID	000327
Solubility (25°C)	6400ng/L
Sensitivity	13.32ng/L
Entrez	406
Expiration Date	Please enquire.