

Product Datasheet

Human Cannabinoid Receptor 1 ELISA Kit (orb1212699)

Description	Human Cannabinoid Receptor 1, CNR1 ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	CNR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,858
Uniprot ID	P21554
Solubility (25°C)	96ng/ml
Sensitivity	0.14ng/ml
Entrez	1268
Expiration Date	Please enquire.