

Product Datasheet

Human Ruvb-like 1 ELISA Kit (orb1212702)

Description	Human Ruvb-like 1, RUVBL1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	RUVBL1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,228
Uniprot ID	Q9Y265
Solubility (25°C)	3200ng/L
Sensitivity	7.15ng/L
Entrez	8607
Expiration Date	Please enquire.