

Product Datasheet

Human Drebrin ELISA Kit (orb1212704)

Description	Human Drebrin, DBN1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	DBN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	71,429
Uniprot ID	Q16643
Solubility (25°C)	3200ng/L
Sensitivity	7.56ng/L
Entrez	1627
Expiration Date	Please enquire.