

Product Datasheet

Human Transcription Factor A ELISA Kit (orb1212716)

Description	Human Transcription Factor A, Mitochondrial, TFAM ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	TFAM
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,097
Uniprot ID	Q00059
Solubility (25°C)	2400ng/L
Sensitivity	5.17ng/L
Entrez	7019
Expiration Date	Please enquire.