

Product Datasheet

Human Tumor Necrosis Factor Ligand Superfamily Member 14 ELISA Kit (orb1212725)

Description	Human Tumor Necrosis Factor Ligand Superfamily Member 14, TNFSF14 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	TNFSF14
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,350
Uniprot ID	O43557
Solubility (25°C)	4800ng/L
Sensitivity	10.92ng/L
Entrez	8740
Expiration Date	Please enquire.