

Product Datasheet

Human Apolipoprotein L1 ELISA Kit (orb1212731)

Description	Human Apolipoprotein L1, APOL1 ELISA Kit
Reactivity	Human
Range	3-900ng/ml
Tested Applications	ELISA
Target	APOL1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,974
Uniprot ID	O14791
Solubility (25°C)	960ng/ml
Sensitivity	1.43ng/ml
Entrez	8542
Expiration Date	Please enquire.