

Product Datasheet

Human Plastin-1 ELISA Kit (orb1212733)

Description	Human Plastin-1, PLS1 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	PLS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	70,253
Uniprot ID	Q14651
Solubility (25°C)	16ng/ml
Sensitivity	0.017ng/ml
Entrez	5357
Expiration Date	Please enquire.