

Product Datasheet

Human Alpha-2C Adrenergic Receptor ELISA Kit (orb1212737)

Description	Human Alpha-2C Adrenergic Receptor, ADRA2C ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	ADRA2C
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,522
Uniprot ID	P18825
Solubility (25°C)	4800ng/L
Sensitivity	12.81ng/L
Entrez	152
Expiration Date	Please enquire.