

Product Datasheet

Human Ephrin B1 ELISA Kit (orb1212746)

Description	Human Ephrin B1, EFNB1 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	EFNB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,007
Uniprot ID	P98172
Solubility (25°C)	32ng/ml
Sensitivity	0.031ng/ml
Entrez	1947
Expiration Date	Please enquire.