

Product Datasheet

Human Campath-1 Antigen ELISA Kit (orb1212753)

Description	Human Campath-1 Antigen, CD52 ELISA Kit
Reactivity	Human
Range	30-9000ng/mL
Tested Applications	ELISA
Target	CD52
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	6,614
Uniprot ID	P31358
Solubility (25°C)	9600ng/mL
Sensitivity	14.96ng/mL
Entrez	1043
Expiration Date	Please enquire.