

Product Datasheet

Human Calcium-dependent Secretion Activator 1 ELISA Kit (orb1212755)

Description	Human Calcium-dependent Secretion Activator 1, CADPS ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CADPS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	152,786
Uniprot ID	Q9ULU8
Solubility (25°C)	3200ng/L
Sensitivity	6.94ng/L
Entrez	8618
Expiration Date	Please enquire.