

Product Datasheet

Human Breast Cancer Anti-Estrogen Resistance Protein 3 ELISA Kit (orb1212758)

Description	Human Breast Cancer Anti-Estrogen Resistance Protein 3, BCAR3 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	BCAR3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	92,566
Uniprot ID	075815
Solubility (25°C)	3200ng/L
Sensitivity	9.33ng/L
Entrez	8412
Expiration Date	Please enquire.