

Product Datasheet

Human Bone Morphogenetic Protein 8A ELISA Kit (orb1212761)

Description	Human Bone Morphogenetic Protein 8A, BMP8A ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	BMP8A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,798
Uniprot ID	Q7Z5Y6
Solubility (25°C)	4800ng/L
Sensitivity	9.02ng/L
Entrez	353500
Expiration Date	Please enquire.