

Product Datasheet

Human Bone Marrow Stromal Antigen 2 ELISA Kit (orb1212762)

Description	Human Bone Marrow Stromal Antigen 2, BST2 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	BST2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,769
Uniprot ID	Q10589
Solubility (25°C)	2400ng/L
Sensitivity	4.72ng/L
Entrez	684
Expiration Date	Please enquire.