

Product Datasheet

Human Beta-Adrenergic Receptor Kinase 1 ELISA Kit (orb1212764)

Description	Human Beta-Adrenergic Receptor Kinase 1, ADRBK1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	GRK2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	79,574
Uniprot ID	P25098
Solubility (25°C)	2400ng/L
Sensitivity	4.95ng/L
Entrez	156
Expiration Date	Please enquire.