

Product Datasheet

Human Neurexin-1-beta ELISA Kit (orb1212771)

Description	Human Neurexin-1-beta, NRXN1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	NRXN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,424
Uniprot ID	P58400
Solubility (25°C)	4800ng/L
Sensitivity	9.91ng/L
Entrez	9378
Expiration Date	Please enquire.