

Product Datasheet

Human Cd83 Antigen ELISA Kit (orb1212776)

Description	Human Cd83 Antigen, CD83 ELISA Kit
Reactivity	Human
Range	1-400ng/L
Tested Applications	ELISA
Target	CD83
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,042
Uniprot ID	Q01151
Solubility (25°C)	480ng/L
Sensitivity	0.58ng/L
Entrez	9308
Expiration Date	Please enquire.