

Product Datasheet

Human Tumor Necrosis Factor Ligand Superfamily Member 13 ELISA Kit (orb1212777)

Description	Human Tumor Necrosis Factor Ligand Superfamily Member 13, TNFSF13 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	TNFSF13
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,433
Uniprot ID	075888
Solubility (25°C)	2400ng/L
Sensitivity	5.61ng/L
Entrez	8741
Expiration Date	Please enquire.