

Product Datasheet

Human Spindlin 1 ELISA Kit (orb1212780)

Description	Human Spindlin 1, SPIN1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	SPIN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,601
Uniprot ID	Q9Y657
Solubility (25°C)	3200ng/L
Sensitivity	8.05ng/L
Entrez	10927
Expiration Date	Please enquire.