

Product Datasheet

Human Complement C1Q Tumor Necrosis Factor-related Protein 4 ELISA Kit (orb1212785)

Description	Human Complement C1Q Tumor Necrosis Factor-related Protein 4, C1QTNF4 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	C1QTNF4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,256
Uniprot ID	Q9BXJ3
Solubility (25°C)	32ng/ml
Sensitivity	0.028ng/ml
Entrez	114900
Expiration Date	Please enquire.