

Product Datasheet

Human B-defensins 1 ELISA Kit (orb1212796)

Description	Human B-defensins 1, DEFB1 ELISA Kit
Reactivity	Human
Range	1-400ng/L
Tested Applications	ELISA
Target	DEFB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	7,420
Uniprot ID	P60022
Solubility (25°C)	480ng/L
Sensitivity	0.58ng/L
Entrez	1672
Expiration Date	Please enquire.