

## Product Datasheet

### Human Zinc Transporter Zip1 ELISA Kit (orb1212797)

|                            |                                                   |
|----------------------------|---------------------------------------------------|
| <b>Description</b>         | Human Zinc Transporter Zip1,<br>SLC39A1 ELISA Kit |
| <b>Reactivity</b>          | Human                                             |
| <b>Range</b>               | 10-2000ng/L                                       |
| <b>Tested Applications</b> | ELISA                                             |
| <b>Target</b>              | SLC39A1                                           |
| <b>Storage</b>             | 2-8°C                                             |
| <b>Note</b>                | For research use only                             |
| <b>Application notes</b>   | Sandwich                                          |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates            |
| <b>MW</b>                  | 34,250                                            |
| <b>Uniprot ID</b>          | <b>Q9NY26</b>                                     |
| <b>Solubility (25°C)</b>   | 2400ng/L                                          |
| <b>Sensitivity</b>         | 5.79ng/L                                          |
| <b>Entrez</b>              | <b>27173</b>                                      |
| <b>Expiration Date</b>     | Please enquire.                                   |