

Product Datasheet

Human Tumor Necrosis Factor-inducible Gene 6 Protein ELISA Kit (orb1212810)

Description	Human Tumor Necrosis Factor-inducible Gene 6 Protein, TNFAIP6 ELISA Kit
Reactivity	Human
Range	0.06-18ng/mL
Tested Applications	ELISA
Target	TNFAIP6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	31,203
Uniprot ID	P98066
Solubility (25°C)	19.2ng/mL
Sensitivity	0.035ng/mL
Entrez	7130
Expiration Date	Please enquire.