

Product Datasheet

Human Tumor Necrosis Factor Superfamily 15 ELISA Kit (orb1212811)

Description	Human Tumor Necrosis Factor Superfamily 15, TNFSF15 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	TNFSF15
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,087
Uniprot ID	O95150
Solubility (25°C)	96ng/mL
Sensitivity	0.155ng/mL
Entrez	9966
Expiration Date	Please enquire.