

Product Datasheet

Human Troponin C ELISA Kit (orb1212818)

Description	Human Troponin C, Skeletal Muscle, TNNC2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	TNNC2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,122
Uniprot ID	P02585
Solubility (25°C)	3200ng/L
Sensitivity	7.17ng/L
Entrez	7125
Expiration Date	Please enquire.